

TECOTHERM NEO

Mobile Therapeutic Hypothermia

- 🕒 Programmed hypothermia protocol
- 🕒 Baby servo controlled
- 🕒 Cools or warms
- 🕒 Compact and lightweight



Total body cooling and warming

Features

TECOTHERM NEO



Three Operational Modes


- Servo-control complete treatment
- Servo-control constant temperature
- Constant mattress temperature



Multiple Mounting Options

Up to 13 mounting options for a custom fit

Small Transport Package



NICU




Air



Ground

Example of Patient Data Collected from Tecotherm Neo



- Infant transfer temperature 35.3°C
- Target temperature 33.5°C (72 hours)
- Target temperature range 33.5°C (+/- 0.5°C)
- Target rewarming temperature 37°C (12 hours/0.29°C per hour)

Order Information

| Part # | Component | Part # | Component |
|----------|---------------------------------|----------|---|
| 731-0154 | Tecotherm Neo with Starter Kit | 731-0147 | Mattress Hose, 2m |
| 731-0255 | Non-Vented Fill-up Set | 731-0148 | Mattress Repair Kit |
| 731-0143 | Interlayer Foils (10pk) | 731-0313 | Reusable Adapter Cable for Skin Temp Probe |
| 731-0144 | Disposable Wrap (3pk) | 731-0152 | Reusable Adapter Cable for Esoph/Rectal Probe |
| 731-0145 | Reusable Wrap | 731-0240 | Disposable Esophageal/Rectal Probes (20pk) |
| 731-0146 | Mattress Hose, 1m | 731-0257 | Reusable Cleaning Loop |
| 731-0258 | Chlorine Dioxide Tablets (30pk) | 731-0312 | Disposable Skin Temp Probe (20pk) |

Europe

International Biomedical AG

Glütschbachstrasse 2, CH-3661 Uetendorf, Switzerland
 T. +41.33.345.66.00 F. +41.33.345.66.08
 welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
 T. 512.873.0033 F. 512.873.9090
 sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.

