





Nasal Cannula | AIRVO 2 | **SPECIFICATIONS**



PRODUCT	OJR416 	OJR418 
Nasal cannula size	L	XL
Cannula weight	13.5 g	13.8 g
<b>COMPATIBLE ACCESSORIES</b>		
F&P Wigglepad™ 2	 WJR112	
<b>MAXIMUM FLOW RATES* (L/min)</b>		
AIRVO 2 	2 - 20	2 - 25
<b>PRODUCT SPECIFICATIONS</b>		
Compatible circuits	900PT531; 900PT561 <sup>1</sup>	
Connection	Easy-click connector	
Compatible humidifier	AIRVO 2 (in Junior mode)	
Quantity	Box of 20	
Box components	F&P Optiflow Junior 2 nasal cannula, user instructions	
<b>PERFORMANCE</b>		
Ambient range °C	18 - 26 °C / 64 - 79 °F	
Humidifier mode	Junior mode	
Usage	Single patient use; Maximum 7 days	
Duration of use	7 days	
Recommended gas source	Medical gas	
<b>COMPONENTS AND COMPOSITIONS</b>		
Predominant materials	Thermoplastic Elastomer; Hydrocolloid; ABS; Stainless steel	
Materials not present	Not manufactured with natural rubber latex, PVC or Phthalates (DEHP, DBP, BBP)	
Manufacturing mode	Non-invasive device; Produced in a controlled working environment	
Disposal	Incineration or according to hospital protocol for cannula; Clamshell packaging and label recyclable PET	
<b>REGULATORY</b>		
Classification	AU-IIa, EU-IIa, Canada-II, US-I	
Country of origin	New Zealand	
Notified body	TÜV SÜD Product Service GmbH CE0123	

\* Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.

<sup>1</sup> Not available for sale in the United States.

Please note that the information in this specifications catalogue (including product information and images) is summarized and provided for illustrative purposes only. Please refer to the relevant user instructions for more information and confirm details with your local Fisher & Paykel Healthcare representative prior to placing an order. Information subject to change without notice.