

# AeroNOx 2.0™

Nitric Oxide Delivery

## Nitric Oxide Delivery System

- ⦿ Meets transport standards
- ⦿ Lightweight and rugged
- ⦿ NO, NO<sub>2</sub> and O<sub>2</sub> monitoring



Where ever you need it, whenever you need it

# Features

# AeroNOx 2.0™

## Nitric Oxide Delivery



### LCD Color Screen

Easy to read and interpret color display.



### Easy Maintenance

Easy access to replace and service the battery and sensors.



### Smart Alarms

Visual and audible alarms updated to the latest alarm standards.

### Transport Mount



NICU



Air



Ground

## Order Information

Part #	Component	Part #	Component
731-0317	AeroNOx 2.0 Nitric Oxide Delivery System	731-9147	INOstat Kit
738-1861	AeroNOx 2.0 Delivery Hose, 3ft, Stainless Steel	731-0274	AeroNOx 2.0 Calibration Kit
888-0115	AeroNOx 2.0 Rechargeable Battery	731-0276	AeroNOx 2.0 Sample/Delivery Kit, Conventional (5pk)
731-0272	AeroNOx 2.0 Transport Mount	731-0277	AeroNOx 2.0 Sample/Delivery Kit, HFV (5pk)
731-9914	AeroNOx 2.0 Pole Clamp Assembly	731-0278	AeroNOx 2.0 Bagger (5pk)

### Europe

#### International Biomedical AG

Glütschbachstrasse 2, CH-3661 Uetendorf, Switzerland  
T. +41.33.345.66.00 F. +41.33.345.66.08  
welcome@int-bio.ch

### Worldwide

#### International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754  
T. 512.873.0033 F. 512.873.9090  
sales@int-bio.com

### int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.

