



TORTLE MIDLINER MAX for NICU and PICU

Provides boundaries and precision head positioning for larger hospitalized infants (Single Patient Use)

MidLiner MAX Benefits:

Supports supine midline /partial lateral position as well as side lying

Velcro front opening

- Provides flexible size adjustment
- Simple and easy to apply which minimizes stress to the baby
- Easy access to the anterior fontanel for ultrasounds and taps
- Easy access for scalp IV's
- Ensures proper fit without slippage
- Decreases intervention time for patients and staff

New design and materials

- Securely holds nasal CPAP interfaces, nasal cannulas, RAM and feeding tubes
- Turtle Aviator Eye shades attach with ease
- MRI, CT and x-ray friendly

Simplifies transition to all positions

- No need to detach or re-apply respiratory or other accessories
- Helps prevent and treat dolichocephaly, plagiocephaly, as well as head preference

Supports "back to sleep" program. Device can be used at all times in the NICU or PICU environment since the patient will always be monitored.



This Turtle is a simple and efficient head repositioning device that will assist staff in providing optimal care for larger infants who need boundaries and optimal head repositioning. Examples could include infants of a diabetic mother, shunted hydrocephalus, NAS infants, and infants with chronic lung disease with long term respiratory care, etc.

Available in three sizes:

SIZE	CASE #	HEAD CIRCUMFERENCE
SMALL	65931	34-42 cms
MEDIUM	65933	42-50 cms

10 units per case.

Includes the following:

- Velcro attachments
- Cannula keys
- Back strap for repositioning



Respiratory Kits:

TYPE	CASE #	COMPATIBLE
Midline	65919	F&P, Draeger, Miniflow, Care Fusion
Goalpost	65921	Hudson, Babi-Plus

Laundry: Can be hand washed and air dried.



CE FDA-cleared